



accesswalk

DESIGNED FOR ALL- TERRAIN MOBILITY

FREEDOM TO MOVE WITH CONFIDENCE

ACCESSWALK Rollator Walker is designed specifically with this in mind. Unlike standard walkers, the ACCESSWALK® is built to navigate challenging surfaces with ease, making it the perfect mobility aid for beaches, parks, or any outdoor adventure.

Benefits of the ACCESSWALK on Sand

Navigating sand can be a challenge for many mobility devices, but the ACCESSWALK Rollator Walker is uniquely suited for this terrain.

Enhanced Traction : The wide, treaded rubber wheels are designed to distribute weight more evenly, preventing the rollator from sinking into soft sand. This ensures smooth movement along beaches or sandy trails.

Stability and Balance: Thanks to its ergonomic design, the ACCESSWALK provides extra stability on uneven surfaces. You won't have to worry about tipping or losing balance, making beach walks enjoyable and stress-free.

Reduced Effort: Unlike traditional walkers that struggle with loose sand, the ACCESSWALK minimizes resistance and glides effortlessly. This reduces the strain on

your arms and legs, allowing for a more comfortable, energy-efficient experience.

Key Features

All-Terrain Wheels: The ACCESSWALK is equipped with large, yellow durable wheels that provide superior traction and stability, even on soft surfaces like sand. Whether you're strolling along the shoreline or exploring a sandy path, our rollator gives you the support and freedom to move confidently.

Lightweight but Sturdy Frame: Crafted with a lightweight yet durable anodized aluminum frame, the ACCESSWALK is easy to maneuver while still providing robust support. This makes it ideal for individuals who need a reliable walker that doesn't feel cumbersome or difficult to transport.

Comfortable Seating and Adjustable Handheld: When you need to rest, the ACCESSWALK has a built-in seat with an adjustable handheld. Whether you're relaxing by the water or taking a break from a beach walk, you'll always have a comfortable place to sit.

User-Friendly Braking System: The intuitive hand brakes ensure that users have complete control, even on slopes or uneven terrain. This makes it easier to stop and start on sandy surfaces, giving you peace of mind while walking.

Foldable and Portable: The rollator easily folds for convenient storage and transport. Whether you're traveling to the beach or heading to a park, the ACCESSWALK is designed to fit effortlessly in your car or other transport.



Specifications

- ▶ **Dimensions of ACCESSWALK shipping in a cardboard box :**
98 cm x 95.5 cm x 28 cm for 18 kg
- ▶ ACCESSWALK fits through standard ADA doorways
- ▶ Lightweight anodized aluminum frame + Front Wheel Stainless Steel casters / axle
- ▶ Foldable for storage and transportation
- ▶ **Folded Up Dimensions :**
(without 12" front wheels & back bar) :
85.1 cm L x 87.5 cm W x 25.4 cm H
- ▶ **Folded Up Dimensions :**
(with 12" front wheels & back bar) :
85.1 cm L x 87.5 cm W x 43.2 cm H
- ▶ **Width of ACCESSWALK :** 87.5 cm
- ▶ **Length :** 76 cm
- ▶ **Net Weight :** 12 kg
- ▶ **Gross Weight :** 18 kg
- ▶ **Maximum weight capacity :** 130 kg
- ▶ Folding seat and textile storage basket
- ▶ **Seat Material :** Polyester mesh incl. foamed cushion
- ▶ **Seat Height :** 59 cm
- ▶ **Seat Dimensions :**
41.5 cm wide x 30 cm deep
- ▶ **Parking brake :**
Yes, activated with the hand brake
- ▶ **Turning Radius Circle :** 90 cm
- ▶ **Yellow "SUNNY" Pneumatic rubber wheels :** Two front wheels of 30.5 cm diameter and two back wheels of 25.4 cm diameter
- ▶ **Adjustable/removable handles from**
82.55 cm to 94 cm high
- ▶ **Width Between Handles :** 45.7 cm
- ▶ **Handles material :** Plastic
- ▶ **Carry Basket :** 5 kg payload
- ▶ **Carry Basket Material :** Mesh fabric
- ▶ **Warranty :** Limited Warranty of 3 Years

www.FUN2ACCESS.com
sales@fun2access.com