

NEED HELP?

**Need Help? Reach out!*



contact@gorampe.com



1-888-957-2673


Name : _____ Address : _____

SLIDING DOOR THRESHOLD - SIDE PERSPECTIVE

CAPTION


- 1** Interior rail
- 2** Exterior Rail
- 6** Free room inside
- 7** Free room outside
- 8** Available width
- 9** Sliding door
- 10** Level

Il manque un espace
 ← → Distance to measure
 - - - Limits between two points



SLIDING DOOR THRESHOLD - SEEN FROM ABOVE

Outside wall



Inside wall

Outside wall



Inside wall



In order to help us provide you with the best product, please send us 2-3 pictures with this form (picture from the inside, outside and general area around the door)

NEED HELP?





Need Help? Reach out!



contact@gorampe.com



1-888-957-2673



MEASUREMENT

3 _____

4 _____

IMPORTANT, measure with a level
 from the highest point of the rail (1-2)

6 _____

7 _____

Make sure there is enough room to get
 on and off the ramp safely

5 _____

Take into account everything that sticks out of
 the wall (mouldings, radiator, etc)

8 _____

Make sure your wheelchair or walker can fit
 through the door
 (TIP: to gain a few inches of clearance, remove
 the stoppers at the extremity of the rail.)